



DONBABY[®]
NURSING SUNDRIES & BABY CARE

PUSH

DOYBABY

Nursing Sundries & Baby Care


Microcomputer Bottle Sterilizer & Dryer

FEATURES

- 1. Only One Touch to sterilize and dry the bottle completely by Microcomputer automatic control.
- 2. Multifunctional pattern control.
- 3. Sterilize 25 minutes to be obtained. High temperature sustained to best efficient bottle sterilizer & dryer.
- 4. The machine is provided with auto-shut off function, simply operation, use at ease during at the same time.
- 5. It could be a single dryer and also to be sterilizing for general accessories for the most products for baby nursing.
- 6. Steam water may make you feel it easy to use it for disinfecting, washing and drying.
- 7. It can take out of the cap of umbrella when in case to reduce the drying function.

COMPONENTS

- 1 Lid handle
- 2 Lid
- 3 Tests tray
- 4 Tongs
- 5 Bottle rack
- 6 Middle base
- 7 Control Button
- 8 Main body base
- 9 Switch



INSTRUCTIONS FOR USE

1. Fill about 200ml-300ml water into stainless steel water tray. (Please read manual carefully before use.)

2. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

3. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

4. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

5. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

6. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

7. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

8. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

9. Press the power button, steam will come down on the bottle rack. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray. When the steam is enough, please take down the bottle rack. When the steam is not enough, please add water into the water tray.

SPECIFICATIONS

Model	MS-13
Product Name	Microcomputer Bottle Sterilizer & Dryer
Power	220V
Power Consumption	150W
Capacity	25 bottles
Material	ABS
Color	Blue
Weight	1.5kg
Dimensions	180x180x180mm
Manufacturer	DOYBABY
Address	1234567890
Phone	1234567890
Fax	1234567890
E-mail	1234567890
Website	1234567890

ISO9001
ISO13485

DOYBABY

Nursing Sundries & Baby Care

Microcomputer Bottle Sterilizer & Dryer

4 procedures

ONE TOUCH

HIGH TEMP. DRYING

AUTO DRY

Auto Shut Off

Free

Disinfecting and drying one touch
25 bottles at a time
Auto shut off

DONBABY

AUTO STEAM STERILIZER

220V/60HZ 60W 1.6A

MADE IN TAIWAN

CE